

Design and Specify your Roam Supply Cart 1.0

Medical Storage Solutions for Healthcare Facilities



solaire
MEDICAL

Roam Supply Cart 1.0 (S3210)

Dimensions: 40.5" wide x 28.75" deep x 80" high

Material: Aluminum

These booklets contain the information you need to design and specify your Roam cart. Once you've completed the design steps, complete the worksheet to summarize your design specifications. Your sales representative and a Solaire Medical Customer Service representative will then complete your order.

For more information about accessories for your Roam cart, see the Accessories Specification Guide under the Resources tab or Accessories under the Products tab on solairemedical.com.



1.0 Cart

The Roam 1.0 cart comes with 1 lockable counter-balanced roll-top door or 2 lockable tempered glass doors. Outfit the interior with shelves, baskets, and trays and the exterior with a variety of side-attached accessories. The cart has 2 interior columns; the column divider can be removed to fit full-width shelves. It has a lightweight aluminum case, durable plastic bumper and base, easy-grip handle, 5" medical-grade casters, and standard key lock.

Follow these steps to design and specify your Roam cart and its interior and exterior accessories.

Step 1

Select interior configuration

Step 2

Select door style

Step 3

Select lock

Step 4

Select surface color


Step 5

Configure interior


Step 6

Configure exterior


Interior configuration



2 18" w compartments



1 36" w compartment



Roam Supply Cart 1.0 (S3210)

Dimensions: 40.5" wide x 28.75" deep x 80" high
Material: Aluminum

Design and Specify your Roam Supply Cart 1.0

Summary of your Roam 1.0 Cart

Interior configuration

- 2 18"w compartments
- 1 36"w compartments

Door style

- Roll-top doors
- Glass doors

Lock

- Keyed lock (standard)
- Electronic lock, roll-top doors
- Electronic lock, hinged doors

Surface color


- White
- Royal
- Red
- Yellow
- Bisque
- Brushed Stainless
- Laminate

Accessorize cart based on configuration and number of compartments

Rows 29 and 30 cannot be used on carts with roll-top doors

30	30
29	29
28	28
27	27
26	26
25	25
24	24
23	23
22	22
21	21
20	20
19	19
18	18
17	17
16	16
15	15
14	14
13	13
12	12
11	11
10	10
9	9
8	8
7	7
6	6
5	5
4	4
3	3
2	2
1	1

Accessorize cart exterior



Left Side

Right Side
















Roam Supply Cart 1.0 (S3210)

Dimensions: 40.5" wide x 28.75" deep x 80" high












Material: Aluminum

Configure interior

Cart with roll-top door has 56 usable cells; cart with glass door has 60 usable cells. Cart without column has 28 usable with roll-top door and 30 usable cells on cart with glass doors

Accessory	Cell use	Quantity
 2" h tray	1 cell	
 4" h tray	2 cells	
 6" h tray	3 cells	
 10" h tray	5 cells	
 6" h basket	3 cells	
 10" h basket	5 cells	
 Standard shelf	1 cell	
 36" w shelf	1 cell	
 Shelf divider	2 cells	
 6" h drawer	4 cells	
 6" h 36" w drawer	4 cells	
 9-hook cath module for 18" w compartment	3 cells	
 9-hook cath module for 36" w compartment	3 cells	
 Pull-out suture model	11 cells	
 Boxed catheter insert	19 cells	

Configure exterior

Accessory	Quantity
 Side accessory rail	
 Glove box	
 Triple glove box	
 Sharps container	
 Flip-up writing surface	
 Waste basket	
 Side storage bin	
 Wire basket	
 Oxygen tank holder	
 Chart holder	
 Embolectomy storage rack	

Note: Sizes are nominal



Storage solutions for healthcare

Solairé Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501
888-435-2256 phone | 616-743-5944 fax | solairemedical.com

©2016 Solairé Medical – A Solairé Portfolio Company, Grand Rapids, Michigan 49501