



Evolve Cabinet with Accessory Pack and Hinged Door

Evolve modular cabinets and casework create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts.

The 84"-high cabinet comes in 3 widths and 2 depths. It has tempered glass doors. Material choices include compression board with a melamine or specified laminate surface, stainless steel, or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. The FlexCell interior provides a durable and cleanable reconfigurable system for trays, baskets, and shelves. Cabinet has a 4"-high toe kick.



Evolve Cabinet with Accessory Pack and Hinged Door



Dimensions (Dimensions are nominal)

Exterior (interior):

18.75"w x 27.25"d x 84"h (14.5"w x 24.75"d x 78.5"h)

26.75"w x 19.25"d x 84"h (22.5"w x 16.75"d x 78.5"h)

36"w x 27.25"d x 84"h (31.75"w x 24.75"d x 78.5"h) (shown)

Cabinet with 6" toe kick is 86"h; cabinet with no toe kick is 80"h

Interior Configuration:

18.75"w cabinet includes 1 single-compartment shelf, 2 shelf dividers, 2 4" trays with dividers, 2 8" trays with dividers, 1 4" basket, 2 8" baskets, and label holders

26.75"w cabinet includes 1 single-compartment shelf, 3 shelf dividers, 2 4" trays with dividers, 2 8" trays with dividers, 1 4" basket, 2 8" baskets, and label holders

36"w cabinet includes 2 single compartment shelves, 4 shelf dividers, 4 4" trays with dividers, 4 8" trays with dividers, 2 4" baskets, 4 8" basket, and label holders

30 FlexCells, 60 in cabinet with center column


Options

- Door can be hinged on right or left side of 18.75"w and 26.75"w cabinet
- Keyed or keyless lock options available
- Optional sloped top available
- Optional 6" toe kick or no toe kick available

Product Number Reference



Product Shown:
SEC2736TGAP



Product Numbers

18.75"w x 27.25"d x 84"h

Compression board: **SEC2719TGLAP** (shown), SEC2719TGRAP

AireCore: SEA2719TGLAP, SEA2719TGRAP

AireCore and Brushed Aluminum: SEB2719TGLAP, SEB2719TGRAP

Stainless steel: SES2719TGLAP, SES2719TGRAP



26.75"w x 19.25"d x 84"h

Compression board: **SEC1927TGLAP** (shown), SEC1927TGRAP

AireCore: SEA1927TGLAP, SEA1927TGRAP

AireCore and Brushed Aluminum: SEB1927TGLAP, SEB1927TGRAP

Stainless steel: SES1927TGLAP, SES1927TGRAP



36"w x 27.25"d x 84"h

Compression board: **SEC2736TGAP** (shown)

AireCore: SEA2736TGAP

AireCore and Brushed Aluminum: SEB2736TGAP

Stainless steel: SES2736TGAP

Evolve Cabinet with Accessory Pack and Hinged Door

Materials and Finishes

Cabinet shell:

Compression board: Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish

Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

Edge banding:

Compression board and AireCore: White ABS

AireCore with Brushed Aluminum finish: Brushed Aluminum acrylic

FlexCell panel:

Injection-molded ABS

Glass-front door:

Compression board: Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim

AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish with tempered glass and clear anodized aluminum trim

Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis and tempered glass

Hinges:

Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Metal concealed European-style

Stainless steel cabinet: Stainless steel

Door pulls:

Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Dull chrome-plated steel bow pull

Stainless steel cabinet: 5" stainless steel bow pull



Storage solutions for healthcare

Solairé Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501
888-435-2256 phone | 616-743-5944 fax | solairemedical.com

© 2018 Solairé Medical, Grand Rapids, MI 49501
© InnerSpace is a registered trademark of Solairé Medical, Grand Rapids, MI 49501
0918