



## Evolve TEE Probe Cabinet with Hinged Door

Evolve modular cabinets and casework create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts.

The 19.25"-deep, 92"-high TEE probe cabinet comes in 2 widths. It has a solid-front or tempered glass door, is vented on the top and bottom, and includes a removable drip tray. Material choices include compression board with a melamine or specified laminate surface, stainless steel, or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. Cabinet has a 4"-high toe kick.



# Evolve TEE Probe Cabinet with Hinged Door



---

## Dimensions (Dimensions are nominal)

Exterior (interior):

**26.75"w x 19.25"d x 92"h (25.25"w x 17"d x 86.5"h)**

**36"w x 19.25"d x 92"h (34.5"w x 17"d x 86.5"h) (shown)**

**92"h cabinet with 6"h toe kick is 94"h**

Interior Configuration:

**26.75"w cabinet holds 3 TEE probes**

**36"w cabinet holds 5 TEE probes**

---


## Options

- Door can be hinged on right or left side of 26.75"w cabinet
- Optional Fan Ventilation Unit available
- Keyed or keyless lock options available
- Optional sloped top available
- Optional 6"h toe kick available

# Product Number Reference



Product Shown:  
**SEC1936PGT**



## Product Numbers

26.75"w x 19.25"d x 92"h

**Compression board:** **SEC1927PGRT** (shown), SEC1927PGLT, SEC1927PSRT, SEC1927PSLT

**AireCore:** SEA1927PGRT, SEA1927PGLT, SEA1927PSRT, SEA1927PSLT

**AireCore and Brushed Aluminum:** SEB1927PGRT, SEB1927PGLT, SEB1927PSRT, SEB1927PSLT

**Stainless steel:** SES1927PGRT, SES1927PGLT, SES1927PSRT, SES1927PSLT



36"w x 19.25"d x 92"h

**Compression board:** **SEC1936PGT** (shown), SEC1936PST

**AireCore:** SEA1936PGT, SEA1936PST

**AireCore and Brushed Aluminum:** SEB1936PGT, SEB1936PST

**Stainless steel:** SES1936PGT, SES1936PST

# Evolve TEE Probe Cabinet with Hinged Door

---

---

## Materials and Finishes

### Cabinet shell:

**Compression board:** Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate

**AireCore:** Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

**AireCore and Brushed Aluminum finish:** Aluminum composite with Brushed Aluminum finish

**Stainless steel:** 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

### Edge banding:

**Compression board and AireCore:** White ABS

**AireCore with Brushed Aluminum finish:** Brushed Aluminum acrylic

### FlexCell panel:

**Injection-molded ABS**

### Solid-front door:

**Compression board:** Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate

**AireCore:** Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

**AireCore and Brushed Aluminum finish:** Aluminum composite with Brushed Aluminum finish

**Stainless steel:** 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis

### Glass-front door:

**Compression board:** Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim

**AireCore:** Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim

**AireCore and Brushed Aluminum finish:** Aluminum composite with Brushed Aluminum finish with tempered glass and clear anodized aluminum trim

**Stainless steel:** 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis and tempered glass

### Hinges:

**Compression board, AireCore, and AireCore with Brushed Aluminum cabinet:** Metal concealed European-style

**Stainless steel cabinet:** Stainless steel

### Door pulls:

**Compression board, AireCore, and AireCore with Brushed Aluminum cabinet:** Dull chrome-plated steel bow pull

**Stainless steel cabinet:** 5" stainless steel bow pull



Storage solutions for healthcare

Solairé Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501  
888-435-2256 phone | 616-743-5944 fax | [solairemedical.com](http://solairemedical.com)

© 2018 Solairé Medical, Grand Rapids, MI 49501  
© InnerSpace is a registered trademark of Solairé Medical, Grand Rapids, MI 49501