

Evolve modular cabinets and casework create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts.



Aluminum channel on glass door

Concealed hinge design

Self-closing door

FlexCell design minimizes space between interior components for maximum storage capacity

FlexCell accepts trays, baskets, and shelves

Seamless edge banding

Aluminum tracks on cabinet interior accept LED lights

Separate toe kick simplifies installation  
Standard 4" h toe kick, 6" h or no toe kick options



Storage solutions for healthcare

Solaire Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501  
888-435-2256 phone | 616-743-5944 fax | [solairemedical.com](http://solairemedical.com)

© 2018 Solaire Medical, Grand Rapids, MI 49501  
© InnerSpace is a registered trademark of Solaire Medical, Grand Rapids, MI 49501  
0918