

Tempo procedure carts are a caregiver's perfect companion. Designed to keep pace with patient care and changing storage needs, Tempo carts are lightweight and easy to maneuver, secure, and adaptable to specific supply and procedure storage.



Accessories attach to side and top rails; can be rearranged easily

Pull-out surface clips into place to prevent sliding when stored within cart

Pull-out surface on 17.75" wide cart extends to left or right, can be easily switched post-purchase

Optional electronic lock design angled to better view and input code

Exterior designed with minimal seams and free of exposed fasteners

Simply attach the rail by sliding it into the side channels. The design prevents rails from being positioned unevenly

Full-width drawer pull lets you grab anywhere along pull to open drawer

Durable one-piece drawer box

Colored drawer inserts add an accent to décor

Glides screw into metal—no slipping or loosening

Durable base takes the bumps and wear of 24/7 use

Casters fit securely in base. Base design allows casters to be easily removed and replaced