


## Tempo 24 Procedure Cart (ST24W5)

Tempo procedure carts are lightweight, easy to maneuver, and adaptable to specific supply and procedure storage.

The Tempo 24 cart has a 3 3", 1 6", and 1 9" drawer configuration. Drawers have soft-return glides and integrated pulls. The seamless plastic work surface has raised edges and an extension that pulls out on either side. The cart accommodates a variety of interior and exterior accessories and is available in six case colors and 12 drawer colors. The cart has a lightweight aluminum case, durable bumper and base, 5" medical-grade casters, and keyed lock.



# Tempo 24 Procedure Cart (ST24W5)



## Materials

Shell: **HDPE core with aluminum exterior**

Base: **16-gauge steel with powder-coat finish**

Corner extrusions: **Anodized aluminum**

Drawers: **18-gauge steel with powder-coat finish**

Work surface: **ABS thermoplastic polymer**

Drawer glides: **Zinc plated, steel ball bearing**

Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

## Dimensions (nominal)

Exterior: **21"d x 29.5"w x 36"h**

Weight: **130 lbs**

## Options

- Electronic keyless lock (SKETC)
- Side-mounted, top-mounted, and interior accessories

